

UK Wheelchairs

Call us
FREE on:
0800 0556377

Free Next Day Delivery in UK only
Mainland Europe Credit Cards Accepted
Prices exclude VAT



£169.99



Just 6.9 KGs
(carry weight)

FEATURES

- Ultra lightweight
- Half folding back
- Quick release rear wheels
- Swing up detachable arm rests
- Flip up detachable footrests
- P U Solid wheelchair tyres
- 12 month warranty
- Optional pneumatic tyres

Free Delivery!

Esteem Eclipse Ultra Lightweight Wheelchair

This superb Eclipse self propelled wheelchair is a very attractive and affordable ultra lightweight solution from Esteem. Thanks to its powder coated aluminium tube frame the Eclipse wheelchair is very light weight making it easy to propel and lift. The uniquely shaped arm rests are designed to flip up with ease to facilitate easy side transfer. The half folding back allows the wheelchair to fit easily in car boots and even on most back seats if necessary. For more detail please visit the web site: www.uk-wheelchairs.co.uk



It has all the standard features including swing away and height adjustable foot rests and puncture proof tyres for peace of mind. Weighing in at 12.7 KGs in total this wheelchair is both lightweight and yet strong. Max user weight 18 Stone/115KG.

Email: sales@uk-wheelchairs.co.uk

UK WHEELCHAIRS

28 Newton Road
Newton Abbot
Devon
TQ12 5AA