

Increase your freedom!

ErgoStand

KP-80



The new “Ergo Stand” is a highly manoeuvrable front wheel drive standing power chair with Karma’s patented “S-Ergo” seating system. Front wheel drive allows the user to negotiate tight spaces indoors. It comes with power backrest recline giving adjustments from 80 to 120 degrees.

The medical benefits of standing

- Activates circulation system and improves cardio-pulmonary function.
- Helps prevent decubitus sores.
- Improves bowel regularity and function.
- Helps bladder and intestinal motility.
- Reduces contractures and ankylosis.

S-Ergo

The turning CURVE of your life



Pressure redistribution





Reduced sliding

Traditional seat



Increased risk of sores due to load on pressure points.



S-Ergo seat



The S-Ergo system (intelligent s-shaped ergonomic seating) provides efficient pressure relief by spreading weight over a greater area, at the same time provides stability and reduces sliding.



100	62	97	44 (40/48)	80°-120°	300Wx2	36Ahx2 50Ahx2	VR 2	62	10°	9km	16km 28km	Steel	135kg	79kg 86kg	55kg

▶ The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 6km/hr (3.7mph). Distance, climbing ability and operating time can vary according to condition. ±1cm ±0.5kg

* KARMA reserves the right to modify information herein without further notice.

Features



Front wheel drive



Power stand



Power recline



Ergo seating



Flip back
Ergonomic
armrests



Knee support



Options



Swing away
joystick



Lights

Sold and Serviced by:

104-40654 Rev.00 SEP.2011

* KARMA reserves the right to modify information herein without further notice. *