



Adjustable footplate height



Flip-back armrest



Quick release rear wheels



Wren 2 (SM-852.2)
▲ Self propel

The new "Wren 2" is a quality lightweight aluminum wheelchair with a lot of adjustments and a very reasonable price.

The "Wren 2" has progressive slowing brakes to allow the attendant to control the wheelchair going down slopes. It also equipped with flip back armrests and detachable swing in/out footrest hangers.

The seat to ground height is adjustable from 19" down to 16.5" on the self propel model, the lowest seat to ground height allows the user to propel the wheelchair with their feet.



▲ Transit

Features

- Aluminium frame
- Seat to ground height adjustable from 19" down to 16.5"
- Flip back armrest
- Detachable swing in/out footrests
- Solid rear wheels (puncture proof)
- Seat sizes 16" x 17" , 17" x 17" , 18" x 17" and 19" x 17"
- Quick release rear wheels (self propel)
- Attendant slowing brakes

Options

- Anti-tips
- Seat belt

Model	items		F/R wheels	Icons representing various wheelchair configurations										Max. User Weight		
	Width															
Wren 2 Q24	16"	41	8"-24"													115kg
	17"	43		59	43	59	32	23	111	48	43	94	14.6kg			
	18"	46		61.5										14.8kg		
	19"	48		64										14.9kg		
Wren 2 F16	16"	41	7"-16"												115kg	
	17"	43		55	43	55	28	23	101	44.5	43	91	13.0kg			
	18"	46		57.5										13.2kg		
	19"	48		60										13.3kg		
																13.4kg

* KARMA reserves the right to modify information herein without further notice *

±1cm ± 0.5kg