



Reno Elite

Captain Seat

Model P325/CS | Power Chair

Reno Elite

Captain Seat

Model P325/CS | Power Chair



Luxury padded Captain Seat.



Large storage basket.



Removable battery pack.



Width adjustable arms.



Easily dismantles without tools.



Rear mud guards.



York Road
Bridgend Industrial Estate
Bridgend
CF31 3TB

Telephone: 01656 674488
Facsimile: 01656 674499
Email: sales@romamedical.co.uk

www.romamedical.co.uk

*Lower limb amputees should refer to the user manual before use.

** Passenger weight, ground conditions and battery condition/size can affect distance of travel.

*** Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.



Special Features

- Luxury padded Captain Seat.
- Easily dismantles without tools.
- Compact and lightweight design.
- Dynamic Shark controller.
- On and off board charging.
- Excellent manoeuvrability.
- Intelligent, regenerative, electromagnetic brakes.
- Large front castors and rear wheels.
- Non marking puncture proof tyres.
- Removable battery pack.
- Rear mud guards.
- Height adjustable, swivel seat.
- Crash tested.

Specifications

Model No.*:	P325
Class (Invalid Carriage):	2
Overall Length:	880mm / 34.6"
Overall Width:	600mm / 23.6"
Max Height (Cptn Seat):	1000mm / 39.3"
Seat Width (Cptn Seat):	440mm / 17.3"
Safe Working Slope:	6 °
Maximum User Weight:	136 Kg / 21.4 stone
Battery Size:	21 Ah
Motor Size:	250 Watt
Range**:	20 Km / 12.5 Miles
Overall Weight (Cptn Seat):	62 Kg / 136.6 lbs
Maximum Speed:	6 Kmph / 4 mph
Wheel Size (Front / Back):	150mm / 250mm
Turning Radius:	550mm / 21.6"
Battery Charger:	2 Amp off board
Ground Clearance***:	50mm / 1.96"
Colours:	Flame Red & Flame Blue

Dealer Stamp