

ULTRA LIGHTWEIGHT ALUMINIUM CONSTRUCTION

WEIGHT
12.2KG

Dash Capri

*Extremely light weight,
efficient opening, easy to store*



www.rhealthcare.co.uk

Rhealthcare

It's ability that matters

Dash Capri

*Extremely light weight, efficient opening,
easy to store*

The Dash Capri with its quality lightweight aluminium frame represents great value to prescribers, carers and users alike. Backed by a 12 month parts guarantee.

Dash Capri sees RHealthcare deliver a fantastic, lightweight, addition to the range of Dash wheelchairs, offering unrivalled value for money combined with high aesthetics and an efficient opening mechanism for ease of lifting and storage.



UP
Brake checks progress



DOWN
Parking brakes engaged



Twin Brake Levers

Makes it easier to apply & release the brake

Ultra light, the Dash Capri folding wheelchair allows the occupant to travel in style. Whether going on holiday, visiting the doctors, out for the day shopping or simply visiting friends, this ultra lightweight (only 12.2kgs) wheelchair is ideal for the occasional user both outdoors and indoors. Complete with efficient opening mechanism and attendant operated brakes.

- Built in tipping lever
- 17" Seat width
- Retarding brake with push to lock brake
- Lightweight aluminium frame (12.2 kgs)
- Padded seat & upholstery for added comfort



The Dash Capri is ideal for users & carers wanting quality, lightness & ease of use



Dash Capri *Technical Information*

WEIGHT INFORMATION

Max occupant weight	100kgs
---------------------	--------

SEATING AREA DIMENSIONS

Width	17"
Depth	16"
Backrest canvas height	460mm
Backrest frame angle from vertical	3 degrees
Seat frame angle from horizontal	3 degrees
Height of armrests	260mm
Footrest to seat adjustable height	370mm - 490mm
Seat front edge to ground	505mm
Push handle height	995mm

CORRIDOR WIDTHS TO TURN (315mm wheel)

Through 360 degrees between walls	1400mm
-----------------------------------	--------

OVERALL DIMENSIONS

Open width	600mm
Open length	1010mm
Open height	995mm
Folded width	260mm
Folded Length	645mm
Folded Height	765mm
Frame wheelbase	400mm

REAR WHEEL

Diameter	300mm
----------	-------

CASTOR

Solid	200mm
-------	-------

As a result of continual product improvement and development, RHealthcare reserve the right to change product specification without prior notice

RHealthcare

It's ability that matters

R Healthcare, Building 2, Philips Campus,
Wellhall Road, Hamilton, ML3 9BZ
Tel: 0845 146 0600 Fax 01246 453 881
sales@rhealthcare.co.uk
www.rhealthcare.co.uk

