



# IM Lightweight Wheelchairs



## i-Explore Wheelchair

The i-Explore is a stylish, robust, lightweight wheelchair which comes with a dual braking system as standard. It also comes with flip back armrests and swing in out footrest hangers. The self propel version has quick release wheels to assist with transportation and storage.

The liquid silver paint finish makes the i-Explore easy on the eye and the aluminium frame together with its thoughtful design makes it easy to use.

### Specifications

MODEL	SIZE INCHES	OVERALL LENGTH CM	OVERALL WIDTH CM	HEIGHT CM	FOLDED SIZE CM	SEAT WIDTH CM	SEAT DEPTH CM	BACK HEIGHT CM
IM-9063	17	105	62	95	30	43	42	46
IM-9063	19	105	68	95	32	48	42	46
IM-9064	17	97	64	95	31	43	42	46
IM-9064	19	97	66	95	33	48	42	46

MODEL	SEAT TO GROUND CM	MAX USER WEIGHT KG	CASTORS SIZE	REAR WHEEL SIZE	G/N WEIGHT KG
IM-9063	49	115	8"	24"	16.0/14.0
IM-9063	49	115	8"	24"	16.0/14.0
IM-9064	49	115	8"	12"	14.5/12.5
IM-9064	49	115	8"	12"	14.5/12.5

Karma Mobility Ltd

Unit 6 Target Park, Redditch, Worcestershire B98 8YN

T: 0845 630 3436 E: info@karmamobility.co.uk

www.karmamobility.co.uk

